







## BDS-MM Family Picosecond Diode Lasers

Optical power up to 60 mW at 50 MHz

- Wavelengths 405, 450, 525, 640, 685, 785, 915 nm
- Power up to 60mW, multi-mode
- Small-size laser module, 40 mm x 40 mm x 120 mm
- Free-beam or multi-mode fibre output
- Pulse repetition rate 20 MHz and 50 MHz, selectable
- Fast on / off / multiplexing capability
- Internal power regulation loop
- All electronics integrated
- No external driver unit
- Simple +12 V power supply
- Compatible with all bh TCSPC devices



Pulse shapes may change due to development in laser diode performance. Power measured in free beam. Coupling efficiency into optical fibres is 60 to 90%, depending on fibre diameter

Designed and manufactured by



Becker & Hickl GmbH  
 Nahmitzer Damm 30  
 12277 Berlin, Berlin  
 Tel. +49 / 30 / 787 56 32  
 Fax. +49 / 30 / 787 57 34  
 email: info@becker-hickl.com  
 www.becker-hickl.com



LASOS Lasertechnik GmbH  
 Carl-Zeiss-Promenade 10  
 07745 Jena, Germany  
 Tel. +49 3641 2944-0  
 Fax +49 3641 2944-17  
 info@lasos.com  
 www.lasos.com



# BDS-MM


## Optical

Repetition Rate, switchable by TTL signal  
 Wavelengths  
 Coupling efficiency into fibres (multi-mode, typical values)  
 Max. optical power  
 Pulse width (FWHM, at medium power)  
 Pulse width (FWHM, at maximum power)  
 Warm-up time for power and pulse shape stabilisation after power on

20 MHz and 50 MHz, other combinations on request  
 405, 445, 525, 640, 685, 785, 915 nm, other on request  
 100µm: 60% 200µm: 80% 500µm: 90%  
 10 to 60 mW at 50 MHz, depends on wavelength version  
 65 to 120 ps  
 120 to 300 ps  
 1 min<sup>1)</sup>


## Trigger Output, to TCSPC Modules

Pulse Amplitude -1V (peak) into 50 Ω  
 Pulse Width 5 ns  
 Leading edge fall time 1 ns, see figure right  
 Output Impedance 50 Ω  
 Connector SMA  
 Jitter between Trigger and Optical Pulse < 10 ps




## Synchronisation Input

Input amplitude +3.3 to +5V into 50 Ω  
 Duty cycle 10 to 30 %. DC equivalent must be < 2.5V  
 Input frequency 20 to 60 MHz  
 Connector SMA  
 Switch between internal clock and sync input automatic, by average voltage at trigger connector



## Control Inputs

Laser ON / Off  
 Response of optical output to on/off signal < 4 us for power 10 to 100%, see figures right  
 External Power Control  
 Response time of optical output to power control < 4 us for power 10 to 100%, see figure right  
 Frequency 50 MHz active H, internal pull-up resistor  
 Frequency 20 MHz active H, internal pull-down resistor  
 Laser runs at 50 MHz with Frequency inputs unconnected



## Power Supply



Power Supply Voltage +9 V to +15 V  
 Power Supply Current at 12V 200 mA to 500 mA 2)

## Mechanical Data

Dimensions 40 mm x 44 mm x 120 mm  
 Mounting holes four holes for M3 screws  
 Heat sink requirements < 2°C / W 3)

## Connector Pin Assignment

Connector version Mini Sub-D  
 Power supply +12V 1, 2  
 GND 4, 5, 9, and case  
 Power control voltage 8  
 Laser On/OFF (active H) 6  
 Frequency 50 MHz (active H, internal pull-up resistor) 7  
 Frequency 20 MHz (active H, internal pull-down resistor) 3



## Maximum Values

Power Supply Voltage 0 V to +15 V  
 Voltage at 'Laser On/Off' and 'Frequency' inputs -2 V to +7 V  
 Voltage at 'Laser Power' input -12 V to +12 V  
 Ambient Temperature 0 °C to 40 °C 3)

- 1) Operation below 13 °C ambient temperature may result in extended warm-up time.
- 2) Depends on case temperature due to laser diode cooling. Cooling current changes with case temperature
- 3) Laser must be mounted on heat sink. Case temperature must remain below 40°C

## Related Products

BDS-SM picosecond diode lasers, BDL-SMN picosecond and CW diode lasers, 375nm, 405nm, 445nm, 473nm, 488nm, 515nm, 640nm, 685nm, 785nm



**Caution: Class 3B laser product. Avoid direct eye exposure. Light emitted by the device may be harmful to the human eye. Please obey laser safety rules when operating the devices. Complies with US federal laser product performance standards.**

## International Sales Representatives



US:  
**Boston Electronics Corp**  
 tcspc@boselec.com  
 www.boselec.com



UK:  
**Photonic Solutions PLC**  
 sales@psplc.com  
 www.psplc.com



Japan:  
**Tokyo Instruments Inc.**  
 sales@tokyoinst.co.jp  
 www.tokyoinst.co.jp



China:  
**DynaSense Photonics Co. Ltd.**  
 info@dyna-sense.com  
 www.dyna-sense.com